

ANKALAPS
EASYCLOSE
Fascia Closure Forceps

**Joystick
Conical Router**




**Fascia
Closure Forceps**



**Double Jaw
Suture Holder**



**Single Jaw
Suture Holder**



 **easy**

 **fast**

 **simple**

 **secure**



Our effort is for human...

www.ankalaps.com

Indications

- ✓ Trocar site closure without extending the incision.
- ✓ Multiple trocar site closure with a single device, per patient.
- ✓ Tacking hernia mesh in position prior to final placement.
- ✓ Ligating abdominal wall bleeders.



Product Groups

Product	Model	Contents
EasyClose	ECS120	Disposable Suture Holder Single Jaw 120 mm
EasyClose	ECS150	Disposable Suture Holder Single Jaw 150 mm
EasyClose	ECS200	Disposable Suture Holder Single Jaw 200 mm
EasyClose	ECD120	Disposable Suture Holder Double Jaw 120 mm
EasyClose	ECD150	Disposable Suture Holder Double Jaw 150 mm
EasyClose	ECD200	Disposable Suture Holder Double Jaw 200 mm
EasyClose System	ECS120S	Disposable Suture Holder Single Jaw 120 mm with Joystick Router
EasyClose System	ECS150S	Disposable Suture Holder Single Jaw 150 mm with Joystick Router
EasyClose System	ECS200S	Disposable Suture Holder Single Jaw 200 mm with Joystick Router
EasyClose System	ECD120S	Disposable Suture Holder Double Jaw 120 mm with Joystick Router
EasyClose System	ECD150S	Disposable Suture Holder Double Jaw 150 mm with Joystick Router
EasyClose System	ECD200S	Disposable Suture Holder Double Jaw 200 mm with Joystick Router